


GluCal In-Vitro Testing System


T. Koeper J. Le Renard Y. Satik Department of Mechanical Engineering | University of California at Santa Barbara

Abstract

Type 1 Diabetes is a chronic illness in which the body cannot regulate its blood glucose levels. Continuous glucose monitors (CGMs) help manage the disease, but they are expensive and painful. Laxmi Therapeutic Devices is developing a new, painless, affordable CGM, but before it can go to human trials it needs to be tested for reliability and accuracy in-vitro. This is where GluCal comes in. We designed GluCal to reproduce the blood glucose profile of a Type 1 Diabetic throughout the It achieves this by mixing two solutions of different concentrations of glucose, to produce any concentration in between. GluCal will be used to validate Laxmi's CGM.

Engineering Specifications	
Max rate of change	5 mg/dL/min
Concentration profile	0-500 mg/dL
No leak	Binary
Flow rate	Run for 24 hours without refill
Fit within a glovebox	1ft x 1ft x 2ft






Optimizing Components

Test Chamber

- Minimizes bubbles by filling from bottom
- Design ensures flow path intersects the sensor without leaks




Pumps

 Miniature peristaltic pumps accurately pump at specified flow rate 10 µL/s


Mixer

 Ensures well-mixed outflow using a stir bar and magnetic stirrer coupled with a DC motor to agitate the solution



Liquid Storage and Waste

1L Pyrex bottles hold sufficient volume for 24-hour experiment



 Fluorescein used as an analog for glucose due to similar diffusion coefficients




Figure 3: Step inputs and measured transient response

significant Absorbance values settle without oscillation, indicating solution is well mixed, and thus Laxmi's CGM can be validated with GluCal.

Acknowledgements

- S. Pennathur, N. Lamb, I. Darke, T. Fork, B. VanAntwerp,
- B. Peterson, T. Marks, T. Susko, S. Laguette, M. Halpin





